

HIV Combo Ag/Ab EIA, 4th Generation Test

West Pacific Medical Laboratory is now performing the HIV Combo Ag/Ab EIA, a 4th Generation HIV assay.

The FDA-approved Fourth Generation GS HIV Combo Ag/Ab EIA is an enzyme immunoassay for the simultaneous qualitative detection of HIV p24 antigen and antibodies to HIV-1 (groups M and O) and HIV type 2 in adult and pediatric populations (children as young as two years old).

Why order this HIV assay:

- ✓ Demonstrated reduction in diagnostic window of 0-20 days, on average¹.
- ✓ Approved by the FDA as a screening method.
- ✓ Intended to assess patient's response to treatments.
- ✓ Early detection can make a positive impact on public health.

A multi-center study conducted to evaluate the GS HIV Combo Ag/Ab EIA, 4th Generation assay concluded that the test significantly reduced the diagnostic window when compared to the 3rd generation screening assay, enabling **earlier diagnosis** of HIV infection¹.

Test Code	Test Name	CPT Codes	Specimen Requirements	Reference Ranges	Turn Around Time
1112	HIV Combo Ag/Ab, Fourth Generation	87389	3mL (2mL) - Serum Separator Tube (SST) Requires SEPARATE Tube	Not Detected	Performed Mon - Fri 24-48 HR TAT

Please contact your WPML Representative for customized test requisition forms.